Swab, Tissue or Fluid?

Unless otherwise specified, multiple tests can not be shared in a single container.

Swab

eSwab

Wound culture with gram stain-LAB503

Anaerobic culture-LAB233 (must be ordered separately)



- Wound culture with gram stain and Anaerobic culture can be performed on a single eSwab
- Do NOT add fluids to eSwab

eSwab

Abscess culture with gram stain (non-superficial infected space)-LAB5000




eSwab

Fungal culture-LAB240



- Send separate swab for Fungal culture



Tissue/Other

Tissue culture with gram stain-LAB271



- Collect in cup with saline

Catheter tip culture-LAB224



- Do not add saline.

Anaerobic culture-LAB233



Fluid

Other Fluids

Aspirate culture with gram stain-LAB5001

- Percutaneous drains
- Needle aspiration of abscess



Urine

IP/ED: Urine culture–AMC/SARAT/GFH LAB5239
UC/OP: Urine culture-LabCorp LAB123301700



Sterile Body Fluid

Body fluid culture with gram stain-LAB269

- Pleural
- Pericardial
- Synovial
- Peritoneal (Ascites)



Respiratory

Respiratory culture with gram stain-LAB900

- Sputum
- Bronchial brushing
- Bronchial washings/BAL



CSF

CSF culture and gram stain-LAB268

- CSF



Anaerobic culture-LAB233

